



Interton Ready

Product Description

The new Interton Ready is a full family of hearing aids built on the GN 2.4 GHz wireless technology featuring ear-to-ear communication, directionality and noise management options for great speech understanding and listening comfort.

The new Interton Ready hearing aids are made for Apple, which allows the user to stream audio directly from Apple devices and the Interton Sound App offers enhanced personal control.


Furthermore, the new Interton Ready family works with the full line of Interton Wireless Accessories for an even better user experience.

The Interton Ready BTE hearing aids are with Nano Coating for optimum durability and meet the IP58 classification for ingress protection.

Models 70 Behind-the-Ear (BTE) hearing aids support closed and open configurations.

The 70 BTE model comes standard with Push Button, Volume Control, Telecoil, and Direct Audio Input (DAI) functionality.

Interton Ready BTE 70 is fully compatible with SureFit™ thin tubes and domes.



© 2018 GN Hearing GmbH. All rights reserved. Interton is a trademark of GN Hearing GmbH. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries.



RD70-DWT

RD70-DW

Model	RD670-DW RD670-DWT	RD470-DW RD470-DWT	RD370-DW RD370-DWT	RD270-DW RD270-DWT
Device Configurations				
Battery size	13			
Colors available	5			
Sound Quality				
"WARP compression (WDRC) - number of channels	14	12	8	6
Comfort				
Adaptive Noise Reduction	3 settings	2 settings	1 setting	1 setting
Adaptive Wind Noise Reduction	2 settings	1 setting	1 setting	-
Environmental Gain Tuner	●	-	-	-
Speech Understanding				
Integrated Directionality	●	-	-	-
Automatic Beamwidth	●	-	-	-
Selectable Beamwidth	●	●	-	-
Auto-Steered Directionality	●	●	●	-
Speech-Focused Directionality	●	●	●	●
Fixed Directionality	●	●	●	●
Feedback Management				
Feedback Manager Plus	●	●	●	●
Auto DFS/Preset Feedback Manager	●	●	●	●
Onboarding				
Acclimatization Manager	●	-	-	-
Convenience				
Ear to Ear Communication (Push Button, Volume Control)	●	●	-	-
AutoPhone	●	●	●	●
Made for Apple	●	●	-	-
Remote Control 2	●	●	●	●
TV Streamer 2	●	●	●	-
Phone Clip 2	●	●	●	-
Multi Mic	●	●	●	-
Micro Mic	●	●	●	-
Interton Sound App/Interton EasyHearing App	●	●	●	-
Robustness				
Nano Coating	●	●	●	●
Fitting Flexibility				
Max Gain Handles	14	12	8	6
Number of Programs	4	4	4	3
Frequency Shifter	●	●	●	●
TSG	●	●	●	●
Interton Fitting 1.0	●	●	●	●
Wireless Fitting	●	●	●	●

400944000GB-18.05-Rev.B

Worldwide Headquarters
Interton A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: +45 4575 1111
interton.com
CVR no. 55082715

United Kingdom
GN Hearing UK Ltd.
Kirtlington Business Centre
Portway, Kirtlington
Oxon OX5 3JA
Tel.: +44 1869 352 800
interton.com

Australia
GN Hearing Australia Pty. Ltd.
Gate C, 19-25 Khartoum Road
Macquarie Technology Park
Macquarie Park NSW 2113
Tel.: (free) 1800 658 955
interton.com

New Zealand
GN Hearing NZ Ltd.
Ground Floor, North Entrance
4 Fred Thomas Drive
Takapuna, Auckland, 0622
Tel.: (free) 0800 900 126
interton.com

Singapore
GN Hearing Pte. Ltd.
2 Kallang Avenue
#07-19 CT HUB
Singapore - 339407
Tel: +65 6320 9388
interton.com

Technical Specifications

		RD70-DWT		
		IEC 60118-0 2nd IEC 711 Ear simulator	IEC 60118-0 3rd IEC 60118-7 ANSI S3.22 2cc coupler	
Reference test gain (60 dB SPL input)	1600 Hz/HFA	45	38	dB
Full-on gain (50 dB SPL input)	Max. 1600 Hz/HFA	62 54	51 48	dB
Maximum output (90 dB SPL input)	Max. 1600 Hz/HFA	131 121	127 116	dB SPL
Total harmonic distortion	500 Hz 800 Hz 1600 Hz	0.5 0.5 0.9	0.2 0.2 0.6	%
Telecoil sensitivity (1 mA/m input)	Max. HFA - SPLIV @ 31.6 mA/m (ANSI)	94		dB SPL
Full-on telecoil sensitivity @ 1mA/m	1600 Hz/HFA	87	80	
Equivalent input noise		25	22	dB SPL
Frequency range (DIN 45605/ANSI)		100-6920	100-6810	Hz
Current drain		1.2	1.2	mA

Data in accordance with IEC60118-0 Edition3.0
2015-06, IEC60118-7 and ANSI S3.22-2009,
supply Voltage 1.3V







Technical Specifications

		RD70-DW		
		IEC 60118-0 2nd IEC 711 Ear simulator	IEC 60118-0 3rd IEC 60118-7 ANSI S3.22 2cc coupler	
Reference test gain (60 dB SPL input)	1600 Hz/HFA	48	43	dB
Full-on gain (50 dB SPL input)	Max. 1600 Hz/HFA	66 58	57 53	dB
Maximum output (90 dB SPL input)	Max. 1600 Hz/HFA	134 126	124 121	dB SPL
Total harmonic distortion	500 Hz 800 Hz 1600 Hz	0.4 1.4 0.9	0.4 0.8 0.7	%
Telecoil sensitivity (1 mA/m input)	Max. HFA - SPLIV @ 31.6 mA/m (ANSI)	98		dB SPL
Full-on telecoil sensitivity @ 1mA/m	1600 Hz/HFA	90	85	
Equivalent input noise		25	20	dB SPL
Frequency range (DIN 45605/ANSI)		100-6560	100-6140	Hz
Current drain		1.2	1.2	mA

Data in accordance with IEC60118-0 Edition3.0 2015-06,
IEC60118-7 and ANSI S3.22-2009, supply Voltage 1.3V







Patents pending

All specifications are subject to change without notice



Patents pending

All specifications are subject to change without notice



Measured according to IEC60118-0 Edition3.0 2015-06 at 1.3 V, impedance 6.2 ohms and 23°C on 2cc coupler. Resp. on 2cc according to IEC60118-7 Second edition 2005-10 and ANSI/ASA S3.22-2009 (HFA average calculated at 1000 Hz, 1600 Hz and 2500 Hz; 0 dB SPL sound pressure equals 20µPa). All measurements without DSP features activated unless indicated otherwise Measurement on O.E.S according to IEC711 1981 According to IEC60118-0 Edition 2 1983 and amendment 1 1994 .

Notes:
O.E.S. = Occluded Ear Simulator
2cc = 2 cm³ coupler
Pi = Acoustic input signal

Basic settings:
Full-on Gain, Reference Test Gain
MPO = Maximum Power Output
Maximum Band Width